

2022



**LAB
PRODUCTS**

Blood Collection Tube



- Ref no:1.12.1
- Specification: 1-10ml
- Material: PET or GLASS
- Packing: 250pcs/bag, 10bag/ctn
- Color option: Transparent
- Different tube available

	Blutentnahmeröhrchen
	Tube de prélèvement de sang
	Tubocolector de sangre
	採血管

Test Tubes



- Ref no:1.12.3
- Material :PS
- Specification:
14*41mm,14*49mm,14*55mm,14*53mm
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Reagenzgläser
	Des tubes à essai
	Tubos de ensayo
	試験管


Pipette



- Ref no:1.12.2.7
- Specification:3ml
- Material :LDPE
- Packing:10000pcs/ctn
- Color option: Transparent

	Pipette
	Pipette
	Pipeta
	ピペット

Microscope Slide



- Ref no:1.12.6
- Specification: 76.2 MM. X 25.4 MM. (3" X 1").
- Thickness: 1 mm. - 1.2 mm.
- Packing: 72pcs/box, 2,880pcs/ctn
- Color option: transparent.

	Objektträger
	Microscope slide
	Portaobjetos de microscopio
	顕微鏡スライド

Urine Tube with Screw Cap



- Ref no:1.12.5
- Specification: 12ml, With a watertight closure
- Material: PS
- Packing:200pcs/bag, 10bags/ctn
- Color option: transparent

	Urinröhrchenmit Schraubverschluss
	Tube d'urine avec bouchon à vis
	Tubo de la orina con tapón de rosca
	スクリューキャップと尿管

Cryovials Tubes



- Ref. no.:1.12.4.14
- 4.5 ml
- Material: polycarbonate
- Specification:conical internal base, o-ring cap
- Packing:500pcs/bag,8bags/ctn
- Color of cap:various for option

	Rohre Kryovials
	Tubes cryovials
	Tubos crioviales
	チューブクライオバイアル

Cryovials Tubes







- Ref. no.:1.12.4.13
- 1.5 ml
- Material: polycarbonate
- Specification: conical internal base, o-ring cap.
- Packing:1000pcs/bag,10bags/ctn
- Color of cap:various for option

	Rohre Kryovials
	Tubes cryovials
	Tubos crioviales
	チューブクライオバイアル

Centrifuge Tube






- Ref. no.:1.12.1.10
- 10ml Centrifuge Tubes
- Conical base
- Packing:250pcs/box, 20boxes/ctn
- Color option: transparent

	Dazu t Rifuge Tsube
	Pour ce faire t Rifuge Tsube
	Tubo de centrifuga
	遠心分離管

PCR Tube



- Ref. no.: 1.12.5
- Size: 0.2ml thin wall
- Material: PS
- Packing:1000pcs/bag,10bags/ctn
- Color option: transparent

	PCR Tsube
	Tsube PCR
	Tubo de PCR
	PCR チューブ

Deep Pore Plate



- Ref. no.: 1.12.16
- Material: PS
- Size : 1.3ml 96-hole deep pore plate
- Specification: Round hole, U-bottom, Sterilization
- Packing: 5pcs/bag, 10bags/ctn
- Color option: transparent

	Porentiefe Platte
	Plaque de pores en profondeur
	Placa profunda de los poros
	ディーブポアプレート

Petri Dish



- Ref no: 1.12.2.1
- Material : High grade PS
- Specification: 60ml, 90ml
- Packing: 1pcs/Plastic bag, 400bags/ctn
- Color option: Transparent

	Petrischale
	Boîte de Pétri
	Placa de Petri
	ペトリ皿

Cuvette Cup



- Ref no: 1.12.4.1
- Material: High grade PS
- Multi-sizes for choice
- Packing: 250pcs/bag, 10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.3
- Material :High grade PS
- Size: Multi-sizes for choice
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.4
- Material :High grade PS
- Specification:TECO Blood Clot Meter
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.5
- Material :High grade PS
- Systems CA530,CA6000,CA100
- Specification: TECO Blood Clot Meter
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.6
- Material :High grade PS
- Specification Systems for Spanish, German Labor
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.7
- Material: High grade PS
- Specification: Trolley for Pacific and German TECO 4
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.8
- Material :High grade PS
- French Organon, zhongshanpeikang
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.9
- Material :High grade PS
- For American Coulter
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.10
- Material :High grade PS
- Specification ForTianjinPulisheng and French STAGO
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.11
- Material :High grade PS
- Specification For beckman
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.12
- Material :High grade PS
- Specification: for Beijing Periong Swiss Serozyme
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.13
- Material :High grade PS
- Specification for Roche
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.14
- Material :High grade PS
- Specification: for Roche, special protein analyzer
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvetasse
	Tazacubeta
	キュベットカップ

Cuvette Cup



- Ref no:1.12.4.15
- Material :High grade PS
- Specification for French Biomerieux
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvelasse
	Tazacubeta
	キュベットカップ

Cuvette Cup






- Ref no:1.12.4.16
- Material :High grade PS
- Specification for American Abbott
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cuvette Tasse
	Cuvelasse
	Tazacubeta
	キュベットカップ

Biochemical Cup





- Ref no:1.12.3.1
- Material :High grade PS
- Specification for Japan Dongya F/CIS series Ø30×70
- Japan Kohden/BotestseriesØ30×60
- Mindray BC-2000 Ø30×52
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Biochemische Cup
	Coupe Biochemical
	Copa bioquímica
	生化学カップ

Biochemical Cup



- Ref no:1.12.3.2
- Material :High grade PS
- Specification for Hitachi 7020,7150 biochemical apparatus
- Beckman CS series
- 700 US Beckman ,Dubang
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Biochemische Cup
	Coupe Biochemical
	Copa bioquímica
	生化学カップ

Biochemical Cup




- Ref no:1.12.3.3
- Material :High grade PS
- Specification: for Hitachi 7020,7150 biochemical apparatus
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Biochemische Cup
	Coupe Biochemical
	Copa bioquímica
	生化学カップ

Biochemical Cup



- Ref no:1.12.3.4
- Material :High grade PP or ABS
- Specification: for swiss rose biochemical apparatus
- For Aulinbas AU Reaction cup
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Biochemische Cup
	Coupe Biochemical
	Copa bioquímica
	生化学カップ

Biochemical Cup



- Ref no:1.12.1.10.1
- Material :High grade PP
- Specification: action cup
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Biochemische Cup
	Coupe Biochemical
	Copa bioquímica
	生化学カップ

Centrifuge Tube



- Ref no:1.12.1.10.3
- Material :High grade PP
- Specification Microcentrifuge tube
- Size:1.5ml
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent (orange,red,yellow,blue)

	Zentrifugenröhrchen
	Centrifugeuse Tube
	Tubocentrífugo
	遠心分離管

Centrifuge Tube



- Ref no:1.12.1.10.4
- Material :High grade PP
- Specification: Microcentrifuge tube
- Size: 0.5ml yellow,1.5ml red ,2ml blue
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Zentrifugenröhrchen
	Centrifugeuse Tube
	Tubocentrífugo
	遠心分離管


Cryovial Tube



- Ref no:1.12.1.12
- Material :High grade PP
- Specification: Cryovial Tube
- Size:5ml blue,1.8ml yellow,1ml red
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Cryovial Rohr
	Tube cryovial
	Tubocriovial
	クライオバイアルチューブ

Racks for Cryovial Tubes



- Ref no:1.12.1.14.1
- Material :High grade PP
- Specification: 36 Weel (5ml)
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Racks fürCryovial Tubes
	Racks fürCryovial Tubes
	Bastidores para criovial Tubos
	クライオバイアルチューブ用ラック

Racks for Cryovial Tubes



- Ref no:1.12.1.14.2
- Material :High grade PP
- 25 position
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Racks fürCryovial Tubes
	Racks fürCryovial Tubes
	Bastidores para criovial Tubos
	クライオバイアルチューブ用ラック

Racks for Cryovial Tube



- Ref no:1.12.114.3
- Material :High grade PP
- Specification: 50 position
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Racks für Cryovial Tubes
	Racks für Cryovial Tubes
	Bastidores para criovial Tubos
	クライオバイアルチューブ用ラック

Racks for Cryovial Tube







- Ref no:1.12.114.4
- Material :High grade PP or PS
- Specification :100 position
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Racks für Cryovial Tubes
	Racks für Cryovial Tubes
	Bastidores para criovial Tubos
	クライオバイアルチューブ用ラック

PCR tube





- Ref no:1.12.1.9.1
- Material :High grade PP
- Specification: PCR Tube 8 strips
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	PCR-Röhrchen
	Tube PCR
	Tubo de PCR
	PCR チューブ

PCR Tube



- Ref no:1.12.1.9.2
- Material :High grade PS
- Specification: PCR Tube 12 strips
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	PCR-Röhrchen
	Tube PCR
	Tubo de PCR
	PCR チューブ

PCR Plate



- Ref no:1.12.1.9.2
- Material :High grade PP
- Specification: PCR plate 96 wells
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	PCR Plate
	Plaque PCR
	Placa PCR
	PCR プレート

Pipette Tips



- Ref no:1.12.13
- Material :High grade PS
- Specification:
5-10ul,200ul,300ul,1000ul,5000ul,
- Material: Medical grade PP
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Pipettenspitzen
	Pointes de pipettes
	Puntas de pipeta
	ピペットチップ

Tip Coupled for Rack



- Ref no:1.12.12.2
- Material :High grade PS
- Specification:
Boxed 200ul tip with filter
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Tip gekoppelt für Rack
	Astuce couplé pour crémaillère
	Tip acoplada para estante
	ヒントラックのために結合されました

Slot for Reagent



- Ref no:1.12.12.5
- Material :High grade PS
- Specification: 25ml
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Slot für Reagenzien
	Emplacement pour réactif
	Ranura para el reactivo
	試薬用スロット

Rack for Tips



- Ref no:1.12.12.2
- Material :High grade PS
- Specification: 1000 up tip box
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Gestell für Tipps
	Rack pour obtenir des conseils
	Rack para consejos
	ヒントについては、ラック

Test Tubes



- Ref no:1.12.3
- Material :High grade PP
- PS material
- Specification: 14*41mm,,14*49mm, 14*55mm, 14*53mm
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

	Reagenzgläser
	Des tubes à essai
	Tubos de ensayo
	試験管

Urine Cup



- Ref no:1.12.1
- Specification:60ml,90ml,120ml
- Materials: cup is made of PP, lid is made of HDPE
- Non-sterile or sterile
- Clear graduation on the cups
- Transparent cup and different color of the lid such as red,
- With label or without label on the cup
- Screw lid
- Packing:250pcs/ctn

	Urinbecher
	Urine tasse
	Taza de la orina
	尿カップ

Stool Container



- Ref no:1.12.2
- Specification:30ml,60ml
- Materials: cup is made of PP, lid is made of HDPE
- Clear graduation on the cups
- With spoon
- Non-sterile or sterile
- Screw lid
- Packing:1000pcs/ctn

	Poo Tasse
	Poo tasse
	Taza de Poo
	プーカップ

Medicine Cup Without Lid



- Ref no: 1.12.3.1
- Specifications: 1oz. (30ml) ,
- Materials: PP
- Non-sterile
- Clear graduation on the cup
- Transparent

	Kleine Medizincup
	Petite tasse de médecine
	Pequeña taza de la medicina
	小さな医学カップ

Medicine Cup Without Lid



- Ref no:1.12.3.2
- Specifications: 2oz.60ml,
- Materials: PP
- Non-sterile
- Clear graduation on the cups
- Transparent

	Kleine Medizincup
	Petite tasse de médecine
	Pequeña taza de la medicina
	小さな医学カップ

Triangular Measuring Cup



- Ref no:1.12.4.1
- PS, material,
- Capacity: 1000ml
- With graduated marks
- Color:transparent

	Triangle Tasse
	Triangle tasse
	la taza del triángulo
	三角形の計量カップ

Triangular Measuring Cup



- Ref no:1.12.4.2
- PP, material,
- Capacity: 1000ml
- With graduated marks
- Color:white

	Triangle Tasse
	Triangle tasse
	la taza del triángulo
	三角形の計量カップ

Test Tube Brush



- Ref no:1.12.5
- Material: Nylon
- Specification: length 20cm,
diameter:2,3,4,5,6,7mm
- Packing: 1pc/bag, 2000pcs/ctn
- Color: white or black

	Reagenzglasbürste
	Brosse pour tubes à essais
	Cepillo para tubos de ensayo
	試験管ブラシ



OFFICE:

**Jumeirah Lake Towers, Cluster G, JBC 1
Office 3801, Dubai, United Arab Emirates**

**info@gsticmedical.com
www.Qcaremd.com**