

2022



# RESPIRATORY PRODUCTS

## Nasal Cannula-Soft Prong



- Ref no:1.14.1.1
- Provide flexibility and a more customized fit for patients
- Color:Clear or green
- Matrial:soft PVC and regular PVC
- Diameter:5mm and 6mm
- Provides crush-proof and kink-proof tubing for continuous flow
- Packing:1pc/polybag,50pcs/ctn or as your required

	Nasenskanüle
	Canule nasale
	Cánula nasal
	鼻カニューレ

## Nasal Cannula-Straight Prong



- Ref no:1.14.1.2
- Color:Clear or green
- Matrial:soft PVC and regular PVC
- Diameter:5mm or 6mm
- Soft tubing to increase patient comfort
- Provides crush-proof and kink-proof tubing for continuous flow
- Packing:1pc/polybag,50pcs/ctn or as your required

	Nasenskanüle
	Canule nasale
	Cánula nasal
	鼻カニューレ

## Nasal Cannula-For Mouth and Nose



- Ref no:1.14.1.3
- Color:Clear or green
- Matrial:soft PVC and regular PVC
- Diameter:5mm or 6mm
- Provide flexibility and a more customized fit for patients
- Soft tubing to increase patient comfort
- Provides crush-proof and kink-proof tubing for continuous flow
- Packing:1pc/polybag,50pcs/ctn or as your required

	Nasenskanüle
	Canule nasale
	Cánula nasal
	鼻カニューレ

## Nasal Cannula-For Mouth and Nose



- Ref no:1.14.1.4
- Color:Clear or green
- Matrial:soft PVC and regular PVC
- Diameter:5mm or 6mm
- Provide flexibility and a more customized fit for patients
- Soft tubing to increase patient comfort
- Provides crush-proof and kink-proof tubing for continuous flow
- Packing:1pc/polybag,50pcs/ctn or as your required

	Nasenskanüle
	Canule nasale
	Cánula nasal
	鼻カニューレ

## Nasal Oxygen Catheter




- Ref no: 1.14.1.5
- Size:FR8,fr10,fr12,fr14
- Suitable for easy and efficient administration of oxygen.
- Soft and smooth open distal end facilitates non-traumatic insertion.
- Funnel connector at proximal end.
- Three type are available,with moveable sponge,with non-moveable
- Sponge at end,and without sponge.

	Nasenskanüle
	Canule nasale
	Cánula nasal
	鼻カニューレ

## O2+CO2 Sampling Cannula



- Ref no: 1.14.1.6
- Adult or Child prong
- Color:Clear or green
- Diameter:5mm or 6mm
- Designed for Oxygen supply and Sampling of exhaled co2 gas simultaneously
- Soft prongs minimize irritation

	Sampling Cannula
	échantillonnage Cannula
	La cánula de muestreo
	サンプリングカニューレ

## Nebulizer Mask



- Ref no: 1.14.2.1
- 6ML capacity anti-spill jar to prevent leakage
- Color: Clear or green
- Diameter: 5mm or 6mm
- Kink-resistant
- Elastic strap fixation
- Latex-free
- Packing: 1pc/polybag, 50pcs/ctn or as your required

	Vernebler Maske
	Masque nébuliseur
	Máscara nebulizador
	ネブライザーマスク

## Nebulizer Mask



- Ref no: 1.14.2.2
- 15ML capacity anti-spill jar to prevent leakage
- Color: Clear or green
- Diameter: 5mm or 6mm
- Kink-resistant
- Elastic strap fixation
- Latex-free
- Packing: 1pc/polybag, 50pcs/ctn or as your required

	Vernebler Maske
	Masque nébuliseur
	Máscara nebulizador
	ネブライザーマスク

## Nebulizer with Mouth Piece



- Ref no: 1.14.2.3
- 15ML or 6 ML capacity anti-spill jar to prevent leakage
- Color: Clear or green
- Diameter: 5mm or 6mm
- Kink-resistant tubing
- Latex-free
- Packing: 1pc/polybag, 50pcs/ctn or as your required

	Vernebler Maske
	Masque nébuliseur
	Máscara nebulizador
	ネブライザーマスク

## Simple Oxygen Masks



- Ref no:1.14.3.1
- Made from non-toxic PVC,clear and soft
- Adjustable nose clip assures comfortable fit
- Adult and Pediatric sized available
- Color: transparent or green
- Diameter:5mm or 6mm
- Packing:1pc/polybag,50pcs/ctn or as your required

	Kontakt xy Gen oder sks
	Masques à oxygène
	Las máscaras de oxígeno
	酸素マスク

## Oxygen Supply Tube



- Ref no:1.14.4.1
- Clear or green color
- Diameter:5mm or 6mm
- Resist kinks and occlusion
- Length can be made as your request
- Packing: 1pc/polybag,50pcs/ctn or as your required

	Kontakt xy Gen Tsube
	Tube d'oxygène
	Tubo de oxígeno
	酸素チューブ

## Non-rebreathing Oxygen Masks



- Ref no:1.14.5.1
- Color:Clear or green
- Diameter:5mm or 6mm
- Safety vent allows for entrainment of room air
- Adjustable nose clip assures comfortable fit
- 500ml-1000ml reservoir bag
- Packing:1pc/polybag,50pcs/ctn or as your required

	Kontakt xy Gen oder sks
	Masques à oxygène
	Las máscaras de oxígeno
	酸素マスク

## Venturi Mask



- Ref no: 1.14.6.1
- Safe, simple delivery of variable oxygen concentrations
- Features 6 color-coded diluters: 24%, 28%, 31%, 35%, 40% and 50%
- Includes adaptor for high-humidity entrainment
- Packing:1pc/polybag,50pcs/ctn or as your required

	Kontakt xy Gen oder sks
	Masques à oxygène
	Las máscaras de oxígeno
	酸素マスク

## Tracheostomy Mask



- Ref no.1.14.7
- For tracheostomy and laryngectomy aerosol therapy. Tubing connector swivels 360 degrees and accepts 22mm I.D. corrugated tubing.
- Soft, clear vinyl mask features a tubing connector for tracheostomy and laryngectomy aerosol therapy.
- Adjustable elastic head strap ensures a comfortable fit.
- Packing:1pc/polybag,50pcs/ctn or as your required

	Tracheotomie Mask
	Trachéotomie Mask
	Máscara de traqueotomía
	気管切開マスク

## Humidifier Bottle



- Ref no:1.14.8
- Disposable, durable plastic with preset audible alarm at 4 -6psi pressure relief
- 250-500 mL usable volume
- Defined maximum and minimum water level lines
- Recessed nipple reduces possibility of breaking off outlet connection
- Packing:1pc/polybag,50pcs/ctn or as your required

	Anfeuchten Flasche
	Humidifier bouteille
	Botella del humidificador
	ウエットボトル

## Oxygen Bag



- Ref no: 1.14.9
- Main specifications : 30L 42L 50L
- That selected by the non-toxic, no harmful chemicals role of the single rubber-coated nylon material from a system of oxygen bags, main hospital for emergency transportation and family beds oxygen
- With clean, high-strength, easy to carry, safe, durable characteristics

	Oxygen bag
	Sac d'oxygène
	Bolsa de oxígeno
	酸素バッグ

## Aerosol Spacer



- Ref no:1.14.10.
- Optimizes the delivery of MDI asthma medication.
- Compatible with most MDI (metered dose inhaler) actuators.
- Helps target medicine to your lungs.
- Helps keep you from spraying drug on the back of your throat.
- Anti Static plastic , Latex free
- With different sizes of mask for pediatric, child, adult use
- Packing: 1pc/box, 50pcs/ctn

	Aerosol Spacer
	Aérosol Spacer
	Aerosol espaciador
	エアゾールスプレーサー

## 3 Ball Breath Exciser



- Ref no:1.14.11
- Wide flow rate range of 600cc to 1200cc per second.
- Single patient use to reduce the risk of cross-contamination.
- Increase your inhaled lung volume and Improve removing mucus or sputum.
- Different color-coded balls in three separate cylinders indicate the approximate flow rate passing through the device.

	Atem exerciser
	Souffle exerciseur
	Ejercitador de respiración
	息エクサ

## PVC Manual Resuscitator



- Ref no:1.14.12
- All the best features at a price that fits your budget.
- Soft cushion mask allows for an effective seal on the patients face.
- The patient port swivels 360 degrees.
- Infant, pediatric and adult sizes available, with masks
- Reservoir style is available in bag or tubing.

	Manuelle Beatmungsbeutel
	Insufflateur manuel
	Reanimación manual
	マニュアル蘇生

## Silicon Manual Resuscitator



- Ref no:1.14.13
- All the best features at a price that fits your budget.
- Soft cushion mask allows for an effective seal on the patients face
- The patient port swivels 360 degrees
- Infant, pediatric and adult sizes available with masks
- Reservoir style is available in bag or tubing.

	Manuelle Beatmungsbeutel
	Insufflateur manuel
	Reanimación manual
	マニュアル蘇生





## **OFFICE:**

**Jumeirah Lake Towers, Cluster G, JBC 1  
Office 3801, Dubai, United Arab Emirates**

**[info@gsticmedical.com](mailto:info@gsticmedical.com)  
[www.Qcaremd.com](http://www.Qcaremd.com)**