

2022



ELECTRODES AND ECG PAPER

Disposable Electrode



- Ref. no. 1.16.1
- Liquid glue or Hardened Hydroge available
- Non-woven Fabrics backing or PE foam available
- Size:Ø55mm
- Form:circular
- Adult
- Packing:50pcs/bag,2000pcs/ctn

	Einweg-Elektrode
	Electrode jetable
	En Dji s Posa b tut Sí filtrado ct
	ディスク電極

Disposable Electrode



- Ref. no. 1.16.2
- Suitable for normal ECG tests in diagnostic or monitoring, it uses an Ag/AgCl sensor element and solid conductive & adhesive Hydro-gel or wet conductive gel for adhesion.
- Liquid glue or Hardened Hydroge available
- Non-woven Fabrics backing or PE foam available

	Einweg-Elektrode
	Electrode jetable
	En Dji s Posa b tut Sí filtrado ct
	ディスク電極

Disposable Electrode



- Ref. no. 1.16.3
- Suitable for normal ECG tests in diagnostic or monitoring, it uses an Ag/AgCl sensor element and solid conductive & adhesive
- Hydro-gel or wet conductive gel for adhesion.
- Silver-Silver chloride, self adhesive, PE foam backing
- Snap fastener, 2.8 x 4.5 cm,
- Packing:50pcs/bag,2000pcs/ctn

	Einweg-Elektrode
	Electrode jetable
	En Dji s Posa b tut Sí filtrado ct
	ディスク電極


Disposable Electrosurgical Pencil



- Ref. no. 1.16.4
- Disposable Electrosurgical Pencil (ESU Pencil) is for general surgery, chest surgery usually.
- High-quality ABS shaped
- 100%% EO gas sterilized packing
- plug standard connector
- 3 meters cable, 2.5m and 5m is for slectable
- Electrode Exposed Length :40mm
- Function: Cut and Coagulation
- Packing:1pc/blister,100pcs/ctn

	Einweg-Elektrode
	Electrode jetable
	En Dji s Posa b tut Sí filtrado ct
	ディスボ電極

ECG Paper



- Ref no:1.16.5
- Different medical paper for EKG /ECG/EEG/CTG/USG checking in Hospital
- Various size available.

	EKG-Papier
	Papier ECG
	Papel del ECG
	ECG 紙



OFFICE:

**Jumeirah Lake Towers, Cluster G, JBC 1
Office 3801, Dubai, United Arab Emirates**

**info@gsticmedical.com
www.Qcaremd.com**