

2022



STERILIZATION POUCH

Heat Sealing Flat Reel



- Ref no. : 1.19.1.1
- Material: 55gsm PE film+60gsm paper
- Size: 50mm*200m ,100mm*200m
can be customized
- Packing: 5-10rolls/ctn, can be customized
- Color option: blue

	Heißsiegel flache Spule
	La chaleur d'étanchéité bobine plate
	carrete de sellado plana de calor
	フラットリールヒートシール

Heat Sealing Gusseted Reel



- Ref no. : 1.19.1.2
- Material: 55gsm PE film+60gsm paper
- Size: 50mm*200m
- Packing: 4rolls/ctn, can be customized
- Color option: blue or green

	Heißsiegel gusseted reel
	Bobine d'étanchéité thermique gousset
	carrete de sellado por calor de fuelle
	ヒートシールマチ付きリール

Self-sealing Flat Pouch



- Ref no. : 1.19.2
- Material: 55gsm PE film+60gsm paper
- Size: 70mm x230mm, 90mm x230mm, 135mm x255mm can be customized
- Packing: 200pcs/inner box, 30 innerboxes/ctn
- Color option: blue

	Selbstdicht flachen Beutel
	Auto-étanchéité poche plate
	Auto-sellado de la bolsa plana
	自己シールフラットポーチ



OFFICE:

**Jumeirah Lake Towers, Cluster G, JBC 1
Office 3801, Dubai, United Arab Emirates**

**info@gsticmedical.com
www.Qcaremd.com**