

2022



SYRINGES AND IV ADMINISTRATION

Disposable Syringe



- Ref no:1.7.1
- Used with hypodermic needles to inject liquid or gases into body tissues
- Luer lock, luer slip
- With needle mounted or by packed, or without needle
- PE/Blister individual package
- The product is EO sterile and has no pyrogen or hemolysis reaction
- 1pc/bag, 100pcs/box or as your required

	Einwegspritze
	Seringue à usage unique
	Jeringuilla desechable
	使い捨て注射器

Disposable Insulin Syringe



- Ref no:1.7.2
- Ultra fine cannula for maximum patient comfort
- Attached needle, no dead space, no waste
- EO STERILIZED
- PE/blister individual package
- Packing: 1pc/bag, 100pcs/box

	Einwegspritze
	Seringue à usage unique
	Jeringuilla desechable
	使い捨て注射器

Irrigation Syringe



- Ref no:1.7.3.1
- Types: Bulb, ring, Flat
- Use for cleaning wounds and getting ready for treatment. Easy way to remove dirt and debris without touching the wound
- The irrigation syringe offers uninterrupted action in a single hand technique for maximum control of Frawn or expelled volume
- The barrel features raised, large and easy to read graduations
- 60ml
- Packing: 1pc/polybag, 500pcs/ctn

	Spülspritze
	Irrigation seringue
	Jeringa de irrigación
	灌漑シリンジ

Disposable Hypodermic Needles





- Ref no:1.7.4
- Used with a syringe to inject substances into the body or extract fluids from it
- High quality stainless steel
- EO-Sterilized,single-use,Non-toxic. Non-pyrogenic
- Sharpness of needle point maximize patient comfort
- 1pc/bag,100pcs/box or as your required

	Injektionsnadeln
	Des aiguilles hypodermiques
	Agujas hipodérmicas
	皮下注射針

Disposable Heparin Cap



- Ref no:1.7.5
- Mainly used with I. V. Catheter together Heparin caps is used to stop blood
- Male luer lock connector
- Latex or polyisoprene
- Color: yellow/transparent
- Packing: blister package, 100Pcs/mid box,5000pcs/ctn
- Carton Meas:45*42*28cm

	Einweg-Heparin Kappe
	Bouchon jetable Héparine
	La heparina casquillo desechable
	使い捨てヘパリンキャップ

Disposable Scalp Vein Set







- Ref no:1.7.6
- Butterfly is connected to the soft non-toxic, non-irritant siliconised tube which does not kink
- The soft wing is easy to grip and fix on patient's body
- Unique protective sheath design for preventing accidental needle sticks
- Color coded wings by size for easy recognition.
- Luer lock and luer slip fitting
- Packing:Blister paper, 2000pcs/carton

	Perfusionsbesteck Set
	Scalp Vein Set
	Conjunto vena del cuero cabelludo
	静脈セットを頭皮

Irrigation Needles



- Ref no:1.7.8
- Multiple angles bended
- Easy for clean out
- Color coded hub by size for clear recognition
- Custom needles can be made to specification
- With/without sharp tip
- Client can choose products according to irrigation part
- Packing:Quantity: 100 pieces/box and 4,000 pieces/carton

	Bewässerung Needles
	Aiguilles d'irrigation
	Agujas de riego
	灌漑ニードルズ

Disposable Blood Collection Needle Holder



- Ref no:1.7.7.2
- Disposable needle holder
- Medical grade PVC
- Latex free;
- Assembly with pen needle or scalp vein set

	Nadelhalter
	Porte-aiguille
	Soporte de la aguja
	持針器

Blood Collection Needle



- Ref no:1.7.7.1
- Color-coded plastic needle hubs
- Provide easy identification gauge of gauges
- Double ended needle may be affixed to major brand needle holders (available separately) with minimum danger of contamination
- Dry, natural rubber sleeve that covers opening on stopper-puncturing end of needle and stops flow of blood as evacuated tube is removed
- Packing:1pc/pe bag,100pcs/box,2000pcs/ctn

	Blutentnahme Nadel
	l'aiguille de prélèvement de sang
	la aguja de extracción de sangre
	採血針

Disposable I.V. Cannula



- Ref no:1.7.11.1
- Specially designed for continuous peripheral intravenous infusion. Transparent flashback chamber allows immediate detection of blood
- Injection port type
- Packing:20pcs/box,720pcs/ctn

	Die intravenöse Katheter
	Un cathéter intraveineux
	Catéter intravenoso
	静脈内カテーテル

Dental Needles



- Ref no:1.7.10
- For Dentist Usage.
- Different color for each model
- Long bevel and short bevel
- Tri-bevel point for easy tissue penetration without coring
- The round cutting edge protecting the gum
- Silicone coating substantially reduces pain and tissue trauma
- Special heat treatment of pipes makes Denject Super-elastic
- Packing:100pcs/box,50box/ctn

	Dentalkanülen
	Aiguilles dentaires
	Agujas dentales
	歯科ニードルズ

A.V. Fistula Needles



- Ref no:1.7.9
- Siliconised ultra thin walled and sharp beveled cannula to minimize trauma to the patient
- Back eye needle to minimize interruption of blood flow
- Flexible butterfly wings for proper fixation
- Color coded wings for easy identification of needle sizes
- Also available with rotating wings
- Tube length: 30 cm
- Sizes available: 15 G, 16 G, 17G
- Packing:1pc/peel bag,100pcs/box, 1000pcs/ctn

	Fistula Nadeln
	Fistule Aiguilles
	Las agujas de la fistula
	瘻ニードルズ

Disposable I.V. Cannula



- Ref no:1.7.11.2
- Specially designed for continuous peripheral intravenous infusion.
- Transparent flashback chamber allows immediate detection of blood
- Pen-holder type
- Packing:20pcs/box,720pcs/ctn

	Die intravenöse Katheter
	Un cathéter intraveineux
	Catéter intravenoso
	静脈内カテーテル

Disposable I.V. Cannula



- Ref no:1.7.11.3
- Specially designed for continuous peripheral intravenous infusion.
- Transparent flashback chamber allows immediate detection of blood
- Butterfly type
- Packing:20pcs/box,720pcs/ctn

	Die intravenöse Katheter
	Un cathéter intraveineux
	Catéter intravenoso
	静脈内カテーテル

Disposable Infusion Set



- Ref no:1.7.12.1
- Spike cap, SYB02 Spike, Small Air inlet (Blue color), #6 Air inlet filter, 4.3x50mm Drip chamber, #12.8 drip chamber filter, Roll regulator, Latex rubber injection (convex protrusion), Luer slip, 21G injection needle, 1.5m tube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Infusionsset
	Perfusion
	Set de infusiones
	輸液セット

Disposable Infusion Set



- Ref no:1.7.12.2
- Spike cap, SYB02 Spike, Small Air inlet (Blue color), #6 Air inlet filter, 14.3x50mm Drip chamber, #12.8 drip chamber filter, Roll regulator, Small Y injection site, Luer slip, 21G injection needle, 1.5m tube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Infusionsset
	Perfusion
	Set de infusiones
	輸液セット

Disposable Infusion Set



- Ref no:1.7.12.3
- Spike with Air inlet and filter, Drip chamber, Roll regulator, Latex rubber injection (non-convex), Luer slip, 22G injection needle, 1.24m tube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Infusionsset
	Perfusion
	Set de infusiones
	輸液セット

Disposable Infusion Set



- Ref no:1.7.12.4
- Spike cap, SYB01 Spike, Small Air inlet (Blue color), #6 Air inlet filter, 14.3x50mm Drip chamber, #13.2 drip chamber filter, Claw Roll regulator, Latex rubber injection (convex protrusion), Luer lock, 21G injection needle, 1.5m tube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Infusionsset
	Perfusion
	Set de infusiones
	輸液セット

Disposable Infusion Set



- Ref no:1.7.12.5
- Spike cap, SYB04 Spike, Small Air inlet (Blue color), #6 Air inlet filter, 15.5x54mm Drip chamber, Roll regulator, Small Y injection site, Luer slip, Luer slip cap, 1.5m elastic and semitransparent tube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Infusionsset
	Perfusion
	Set de infusiones
	輸液セット

Disposable Infusion Set



- Ref no:1.7.12.6
- Spike cap, SYA03 Spike,14.3x50mm Drip chamber, #12.8 drip chamber filter, Aluminum clip, Small Y injection site, Luer slip, Luer slip cap, 1.5m tube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Infusionsset
	Perfusion
	Set de infusiones
	輸液セット

Disposable Infusion Set



- Ref no:1.7.12.7
- Spike cap, SYB06 Spike, Big Air inlet (Blue color), #8 Air inlet filter, 6.3x64mm Drip chamber, #14.2 drip chamber filter, Claw Roll regulator, Small Y injection site, Luer slip, 21G injection needle, 1.5m tube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Infusionsset
	Perfusion
	Set de infusiones
	輸液セット

Infusion Set with Burette



- Ref no:1.7.13
- Spec: 100ml, 150ml (Hard Type)
- Vented or ABS(White or Milk-white)
- Top Flow Regulator--PE
- 100ml or 150ml Graduated Burette, with Float Valve inside
- Drip Chamber--PVC, Transparent, with Medication Filter
- 60 drops = 1±0.1ml
- Tubing--PVC,Transparent
- Screw Adapter--ABS,Transparent
- Packing:1pc/bag,50pcs/box or as your required

	infusionsset
	Perfusion
	Set de infusiones
	輸液セット

Continu-FLO Solution Set



- Ref no: 1.7.15
- Specification: luer lock, one spike, one champer,Two injection sites, one slide clamp, tube 2.6m,10 drops/ml, roller clamp with air vent
- Packing: 1pc/bag, 20pcs/big bag, 400pcs/ctn

	Continu-FLO Lösung Set
	Continu-FLO Set
	Continu-FLO conjunto de soluciones
	Continu-FLO ソリューションセット

Secondary Line







- Ref no: 1.7.16
- Spike,
- One champer,
- Transparent tube 1.5m,
- 10drops/ml,
- Roller clamp with air vent
- Without extension hanger
- Packing:1pc/bag, 25pcs/big bag, 400pcs/ctn

	Nebenleitung
	ligne secondaire
	línea secundaria
	二次ライン

Blood Transfusion Set




- Ref no:1.7.14.1
- Spike cap, SXB01 Spike, Big Air inlet, #8 Air inlet filter, 17.3x98mm Drip chamber, 17.3x62mm Drip chamber filter, Roll regulator, Latex rubber injection (convex protrusion), Luer slip, 18Gx1 " needle, 1.5m tube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Blutübertragung
	Transfusion sanguine
	Transfusión de sangre
	輸血セット

Blood Transfusion Set



- Ref no:1.7.14.2
- Spike cap, SXA01 Spike, 17.3x98mm Drip chamber, 17.3x62mm
- Drip chamber filter, Roll regulator, Latex rubber injection (convex protrusion), Luer slip, 18Gx1 " needle,1.5mtube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Blutübertragung
	Transfusion sanguine
	Transfusión de sangre
	輸血セット

Blood Transfusion Set






- Ref no:1.7.14.3
- Spike cap, SXB01 Spike, Big Air inlet, #8 Air inlet filter, 17.3x98mm Drip chamber, 17.3x62mm Drip chamber filter, Roll regulator, Big Y injection site, Small Luer lock, 18Gx1 " needle, 1.5m tube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Blutübertragung
	Transfusion sanguine
	Transfusión de sangre
	輸血セット

Blood Transfusion Set



- Ref no:1.7.14.4
- Spike cap, SXA01 Spike, 17.3x98mm Drip chamber, 17.3x62mm Drip chamber filter, Claw Roll regulator, Small Luer lock, 18Gx1 " needle, 1.5m tube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Blutübertragung
	Transfusion sanguine
	Transfusión de sangre
	輸血セット

Blood transfusion set



- Ref no:1.7.14.5
- # 20 Spike needle, Soft Drip chamber, Drip chamber filter, Claw Roll regulator, Latex rubber injection (convex protrusion) Luer slip, 18Gx1 " needle, 1.5m tube
- Packing:1pc/polybag,20pcs/mid bag,500pcs/ctn

	Blutübertragung
	Transfusion sanguine
	Transfusión de sangre
	輸血セット

Three-way Stopcock



- Ref no:1.7.17
- Minimum priming volume required for accurate drug administration
- Rotating Male Luer and two Female Luer connectors
- Luer lock/Luer slip
- White/red/blue
- With extension tube,the tube length according to customerequirements
- Packing:1pc/bag,50pcs/box or as your requiremen

	Absperrhahn
	Robinet d'arrêt
	Llave de cierre
	栓

Extension Pressure Line



- Ref no:1.7.18
- Considerably reduces the chances of mechanical irritation and infection thereby reducing trauma to the patient.
- Extends the pathway between the infusion site and the source of infusion
- Packing:1pc/bag,100pcs/box or as your required

	Verlängerungsleitung
	ligne d'extension
	línea de extensión
	内線

Single Blood Bag



- Ref.no. 1.7.19.1
- Rolled single blood bag
- Material: PVC, medical grade
- Size:250ml,350ml,450ml,500ml
- Needle Gauge:16G or 17G
- Packing:100pcs/ctn

	Blutkonserve
	Sac de sang
	Bolsa de sangre
	血液バッグ

Double Blood Bag



- Ref.no. 1.7.19.2
- Rolled double blood bag
- Material: PVC, medical grade
- Size:250ml,350ml,450ml,500ml
- Needle Gauge:16G or 17G
- Packing:100pcs/ctn

	Blutkonserve
	Sac de sang
	Bolsa de sangre
	血液バッグ

Triple Blood Bag



- Ref.no. 1.7.19.3
- Rolled triple blood bag
- Material: PVC, medical grade
- Size:250ml,350ml,450ml,500ml
- Needle Gauge:16G or 17G
- Packing:100pcs/ctn

	Blutkonserve
	Sac de sang
	Bolsa de sangre
	血液バッグ

Quadruplet Blood Bag



- Ref.no. 1.7.19.4
- Material: PVC, medical grade
- Size:250ml,350ml,450ml,500ml
- Needle Gauge:16G or 17G
- Packing:100pcs/ctn

	Blutkonserve
	Sac de sang
	Bolsa de sangre
	血液バッグ

Central Venous Catheter Kit



- Ref no:1.7.20
- Excellent soft-tip to significantly reduce risk of vascular trauma.
- Medical grade material: catheter automatically softens under the temperature of 37°.
- With radiopaque materials, ensure correct placement of the catheter tip.
- 1set/blister,20sets/ctn

	Verlängerungsleitung
	ligne d'extension
	línea de extensión
	内線

Needle For Anesthesia



- Ref. no.1.7.21
- Full sizes of anesthesia needles
- Identified as spinal needle bevel quincke tip, pencil point tip and epidural needle
- 304 stainless steel needle tubes and stylets
- Packing:1pc/peel bag,1000pcs/ctn

	Anesthesia Nadel
	Aiguille d'anesthésie
	Hermanas Teshia Hey d tut
	麻醉針

Combined Spinal and Epidural Kit



- Ref. no.1.7.22.1
- Epidural needle(tuohy needle), 16/18G*80mm
- Spinal needle:25G*110mm,pencil point
- Epidural catheter:multiport with mark 0.8/1.0mm≥850mm,
- Epidural filter:0.22um,hydrophobic membrane.Luer lock
- LOR syringe,7ml,10ml available.
- Catheter adaptor:TBA,secure the catheter completely.
- Packing:10set/box,100set/ctn

	Spinal und epidurale Kit
	Spinal et Kit périurale
	S Pina l Anne d Epidzura l-t a cabo
	脊髄と硬膜外キット

Epidural Kit



- Ref. no.1.7.22.2
- Epidural puncture needle 16G/18G
- Epidural anaesthesia tube \varnothing 1/0.8mm
- LOR syringe 7ml/10ml
- Tube connector \varnothing 1.0mm
- Epidural filter 0.22 micron
- Packing: 10set/box,100set/ctn

	Epidurale Kit
	Epidzura l-t out
	Epidzura l-t a cabo
	硬膜外キット

Spinal Kit



- Ref. no.1.7.22.3
- Spinal needle:25G*90mm,
- Epidural catheter:multiport with mark 0.8/1.0mm≥850mm,
- Epidural filter:0.22um,hydrophobic membrane.Luer lock
- LOR syringe,7ml,10ml available.
- Catheter adaptor:TBA,secure the catheter completely.
- Packing: 10set/box,100set/ctn

	Rücken Ki
	Spinale Ki
	Kit de médula
	脊髄キット



OFFICE:

**Jumeirah Lake Towers, Cluster G, JBC 1
Office 3801, Dubai, United Arab Emirates**

**info@gsticmedical.com
www.Qcaremd.com**