

2022



**PLASTIC  
DISPOSABLE**

## Medical Bowls



- Ref no:1.8.1.1
- Specifications: 500ml,1000ml
- Materials: PP
- Transparent
- With clear graduation on
- Packing:500pcs/ctn

	Medizinische Handgelenk
	Petite tasse de médecine
	Muñeca médica
	医療手首

## Medical Bowls



- Ref no: 1.8.1.2
- Specifications: 120ml, 250ml, 500ml, 1500ml,8500ml
- Materials: PP
- Green
- With clear graduation on
- Packing can be customized

	Medizinische Handgelenk
	Petite tasse de médecine
	Muñeca médica
	医療手首

## Square Basin



- Ref no:1.8.1.3
- PP material, L32x W25x H11cm
- Durable and strong
- Color can be customized
- Packing can be customized

	Quadratische Becken
	Bassin carré
	Cuenca del cuadrado
	方盆

## 200ml Kidney Tray



- Ref no:1.8.1.4.1
- PP material
- Size: L16x W8x H3.2cm 200ml
- Ideal for contain fluid, measuring output, routing oral care -With graduated marks
- Color: green
- Packing :250pcs/ctn

	Nierenschale
	Haricot réniforme
	Batea
	腎臓トレイ

## Kidney Tray







- Ref no:1.8.1.4.2
- PP material
- Size:500ml 700ml
- Ideal for contain fluid, measuring output, routing oral care
- With graduated marks
- Color: blue, rose, yellow, green
- Packing:250pcs/ctn

	Nierenschale
	Haricot réniforme
	Batea
	腎臓トレイ

## Urinal for Men with Cap







- Ref no:1.8.1.5.1
- Size:12.5\*9\*25cm
- PE material,
- With Lid
- Capacity: 1000cc
- Packing :50pcs/ctn

	Männliche Urinflasche
	Bouteille d'urine mâle
	Botella de orina masculina
	男性の尿ボトル

## Urine Bottle for Male







- Ref no:1.8.1.5.2
- Weight:125g/pc
- PP material,
- With Lid
- Capacity: 1000cc
- Packing:100pcs/ctn

	Männliche Urinflasche
	Bouteille d'urine mâle
	Botella de orina masculina
	男性の尿ボトル

## Urine Bottle for Male







- Ref no:1.8.1.5.3
- Weight:155g/pc
- PC material,
- With Lid
- Capacity: 1000cc
- Packing:100pcs/ctn

	Männliche Urinflasche
	Bouteille d'urine mâle
	Botella de orina masculina
	男性の尿ボトル

## Urinal for Female



- Ref no:1.8.1.6
- Size:12.5\*9\*25cm
- PE material,
- With Lid
- Capacity: 1000cc
- Packing:100pcs/ctn

	Weibliche Urinflasche
	Bouteille d'urine féminine
	Botella orina femenina
	女性の尿ボトル

## Container Pan



- Ref no:1.8.1.7.1
- PP material,16.5\*11.5\*9.5cm
- With double graduated marks
- Capacity:900cc
- Latex free
- Packing:100pcs/ctn

	Containerwanne
	Pan de récipient
	Pan de contenedores
	女性の尿ボトル

## Container Pan



- Ref no:1.8.1.7.2
- PE material,16.5\*11.5\*9.5cm
- With double graduated marks
- Capacity:1000cc
- Latex free
- Packing:100pcs/ctn

	Containerwanne
	Pan de récipient
	Pan de contenedores
	女性の尿ボトル

## Fracture Bedpan







- Ref no:1.8.1.9.1
- PP material
- Size:425\*275\*104mm
- Color:Gray ,Gold
- Packing:12pcs/ctn

	Bettpfanne
	Bassin
	Calentador de cama
	便器

## Fracture Bedpan



- Ref no:1.8.1.9.2
- PP material
- Size:300\*230\*8mm
- Color:Gray ,Gold
- Packing:12pcs/ctn

	Bettpfanne
	Bassin
	Calentador de cama
	便器

## Stackable Bed Pan



- Ref no:1.8.1.10
- PP material
- Size:38\*31\*12cm
- Stackable
- Color: rose, yellow, blue
- Packing:12pcs/ctn

	Bettpfanne
	Bassin
	Calentador de cama
	便器

## Disposable Plastic Pitcher



- Ref no:1.8.1.11
- PP material,
- With Lid
- Capacity:1000cc
- Packing:50pcs/ctn

	Kunststoff-Krug
	Pitcher plastique
	La jarra de plástico
	プラスチックピッチャー

## Drinking Cup



- Ref no:1.8.112
- Specifications: 8oz. 240ml,
- Materials: PP OR PS
- Non-sterile
- Clear graduation on the cups
- Packing can be customized

	Kleine Medizincup
	Petite tasse de médecine
	Pequeña taza de la medicina
	小さな医学カップ

## Bedside Tray



- Ref no:1.8.113.1
- PP material,
- Size: 213\*137mm
- Color: Graphite,pink
- Packing:500pcs/ctn

	Nacht Tray
	Plateau de chevet
	Bandeja de cabecera
	ベッドサイドトレイ

## Bedside Tray Wiridges



- Ref no:1.8.113.2
- PP material,
- Size: 250\*150mm
- Color: Graphite,pink
- Packing:100pcs/ctn

	Nacht Tray
	Plateau de chevet
	Bandeja de cabecera
	ベッドサイドトレイ

## Beside Tray



- Ref no:1.8.113.3
- PP material,
- Size: 200 x 75 x30mm
- Weight: 20g/pc
- Color: green
- Packing can be customized

	Nacht Tray
	Plateau de chevet
	Bandeja de cabecera
	ベッドサイドトレイ

## Plastic Tray



- Ref no:1.8.113.4
- Material: PP
- For single only use
- Type:3 compartment tray
- Packing can be customized

	Nacht Tray
	Plateau de chevet
	Bandeja de cabecera
	ベッドサイドトレイ

## Plastic Tray



- Ref no:1.8.113.5
- Material: PP
- For single only use
- Type:full deep tray
- Packing can be customized

	Nacht Tray
	Plateau de chevet
	Bandeja de cabecera
	ベッドサイドトレイ



## Plastic Tray



- Ref no:1.8.113.6
- Material: PP
- For single only use
- Type:halfdeep tray
- Packing can be customized

	Nacht Tray
	Plateau de chevet
	Bandeja de cabecera
	ベッドサイドトレイ

## Plastic Tray



- Ref no:1.8.113.7
- Material: PP
- For single only use
- Type: 2 compartment tray
- Packing can be customized

	Nacht Tray
	Plateau de chevet
	Bandeja de cabecera
	ベッドサイドトレイ

## Plastic Tweezer Forcep



- Ref no:1.8.114.1
- PP material
- Non-toxic, Odorless
- Size:12.8cm
- Color:green or blue
- Packing:1000pcs/ctn

	Kunststoff forcep
	Forcep plastique
	Pinza de plástico
	プラスチック鉗子


## Plastic Tweezer Forcep



- Ref no:1.8.114.2
- PP material
- Non-toxic, Odorless
- Size:10cm
- Color: black
- Packing:2000pcs/ctn

	Kunststoff forcep
	Forcep plastique
	Pinza de plástico
	プラスチック鉗子

## Plastic Tweezer Forcep



- Ref no:1.8.114.3
- PP material
- Non-toxic, Odorless
- Size:12.7cm
- Color: blue
- Packing:1500pcs/ctn

	Kunststoff forcep
	Forcep plastique
	Pinza de plástico
	プラスチック鉗子

## Plastic Tweezer Forcep



- Ref no:1.8.114.4
- PP material
- Non-toxic, Odorless
- Size:13 cm
- Color: blue
- Packing :1000pcs/ctn

	Kunststoff forcep
	Forcep plastique
	Pinza de plástico
	プラスチック鉗子

## Plastic Dressing Forcep




- Item no :1.8.114.5
- Material: ABS, Nylon
- Non-toxic, Odorless
- Size:19cm,24cm
- Packing:1000pcs/ctn

	Kunststoff forcep
	Forcep plastique
	Pinza de plástico
	プラスチック鉗子

## Soap Dish



- Item no :1.8.115.1
- Material:PP
- Size: 3.75X2.75X 0.5Inch
- Color : Graphite, White. Pink
- Packing:100pcs/ctn

	Seifenschale
	Antenne SOAP
	Jabonera
	ソープディッシュ

## Soap Dish







- Item no : 1.8.115.2
- Material:PP
- Size : 12.5\*7cm,
- Color :blue
- Packing:500pcs/ctn

	Seifenschale
	Antenne SOAP
	Jabonera
	ソープディッシュ

## Gail Pot




- Item no : 1.8.116
- Material:PP
- Capacity:60ml
- Color: transparent
- Packing can be customized

	Flaschen Medizin
	Bouteilles de médecine
	Botellas de la medicina
	医学のボトル

## Amber Bottle



- Item no : 1.8.117
- Body material:PS Cap material:HDPE
- Capacity:30ml
- Color: amber
- Packing:500pcs/ctn

	Amber-Flasche
	Bouteille ambre
	Frasco de color ámbar
	アンバーボトル

## Plastic Spoon



- Item no : 1.8.118
- Material:PP
- Different size available
- Different color available
- Packing can be customized

	Löffel
	Cuiller
	Cucharas
	スプーン

## Plastic Knife



- Item no : 1.8.1.19
- Material:PP
- Different size available
- Different color available
- Packing:1000pcs/ctn

	Kunststoff-Messer
	Couteau en plastique
	Cuchillo de plástico
	プラスチック製のナイフ

## Plastic Fork




- Item no : 1.8.1.20
- Material:PP
- Different size available
- Different color available
- Packing can be customized

	Plastikgabel
	Fourchette en plastique
	Tenedor de plástico
	プラスチックのフォーク

## Ear Syringe



- Item no : 1.8.1.21
- Material:PVC
- Size:30ml,60ml
- Different color available
- Packing can be customized

	Ear Waschball
	Ear lavage balle
	Oreja de la bola de lavado
	耳ウォッシュボール

## Probe Covers







- Ref. no.: 1.8.1.22
- Ear thermometer probe covers
- Use only with Thermoscan pro series, instant thermometers
- Packing: 20pcs/box, 1000box/ctn

	Sondenüberzüge
	Couvercles de sonde
	Cubiertas de muestra
	プローブカバー

## Medical Dispenser



- Ref no: 1.8.1.23
- Size: 22\*6.5\*1.4cm
- Materials: PP
- Non-sterile
- Color: Blue, White

	Medizinische Spender
	Appareil de distribution médical
	Dispensador médica
	医療ディスペンサー

## Pill Box



- Ref no: 1.8.1.24
- Size: 16.5\*11\*2.2cm
- Materials: PP
- With 28 Lattice

	Pill-Box
	Boîte à pilules
	Caja de pastillas
	印籠

## Pill Box



- Ref no: 1.8.1.25
- Size: 20.9\*12.3\*2.1cm
- Materials: PP
- With 28 Lattice

	Pill-Box
	Boîte à pilules
	Caja de pastillas
	印籠

## Sanitizer Dispenser



- Ref no: 1.8.1.26
- Capacity: 1000ml
- Color: White
- Powder: 1.2W
- With Stand or without stand

	Sanitizer Dispenser
	Sanitizer Dispenser
	Dispensador del desinfectante
	消毒ディスペンサー

## Plastic Towel Clamp



- Ref no: 1.8.1.27
- Material: ABS
- Color: White

	Kunststoff-Tuch-Clamp
	Pince-serviettes en plastique
	Toalla abrazadera plástica
	プラスチック製タオルクランプ



## **OFFICE:**

**Jumeirah Lake Towers, Cluster G, JBC 1  
Office 3801, Dubai, United Arab Emirates**

**[info@gsticmedical.com](mailto:info@gsticmedical.com)  
[www.Qcaremd.com](http://www.Qcaremd.com)**