

2022




SURGICAL DRAPES AND PACKS

Two Layer Folio Drape (Biplex Drape)



- Ref no.: 1.9.1.1
- Material: PE film+ viscose or PP+PE or customized
- Size: 50*50cm,45*75cm or customized.
- Packing: customized
- Style: disposable, sterile
- Color option: blue

	Zweischichtige foliant Drapierung
	Deux couches drapé folio
	Para Aie r D Foio la Paix
	層フォリオドレーブ


Field Drape



- Ref no.: 1.9.1.2
- Material:16gsm paper+10gsm PE+16gsm paper.
- With or without hole
- Use hot-pressing technology.
- Size: 45*65cm or customized
- Packing: customized
- Style: disposable ,sterile

	Feld drapieren
	Champ drapé
	Campo drapeado
	フィールドドレーブ


Paper+PE



- Ref no.: 1.9.1.3
- Material: Paper+PE or 2 layer Paper+PE
- With or without hole (circle, Rhombus, Rectangle)
- With or without adhesive tape
- Color: white or blue
- Style: disposable, sterile

	Papierbeschichtung
	Paper Coating
	Recubrimiento de papel
	紙加工


Three Layer Drape Sheet (Triplex Drape)



- Ref no.:1.9.1.4
- Material: Hydrophilic PP+ meltblown fabric +PE film or customized
- Size: 35*50cm,50*50cm or customized
- Packing: customized
- Style: disposable, sterile
- Color option: blue

	Dreischichtige Drapierung Blatt
	Trois couches Drapé Sheet
	Th Rey Aie r D la Paix Sheet
	三層ドレープシート


Arthroscopy Drapes (Knee Drape)



- Ref no.: 1.9.1.5
- Material: 45gsm SMS with liquid collection pouch
- FDA Approved
- Size : 200*300cm or customized
- Packing: 20pcs/ctn
- Style: disposable, sterile
- Color option: blue

	Arthroskopie Drapierung (Knie Drapierung)
	Trou de serviettes de l'articulation du genou
	Toalla de la articulación de la rodilla
	ひざ関節洞

Ophthalmic Drape



- Ref no.:1.9.1.6
- Material:40gsm SMS
- By CE/ISO Approved
- Packing: 20pcs/ctn
- Style:disposable ,sterile
- Color option: blue

	Augenheilkunde. Drapierung
	Trou de tissu ophtalmique
	Toalla de oftalmología
	眼科洞

Split Drape (U Drape)



- Ref no.:1.9.1.7
- Material:28gsm viscose+ 27gsm PE film
- Style: disposable, sterile
- High-protective
- Size : 80*100CM
- Packing: 1PC/POUCH, 100PCS/CTN

	U - Drapierung
	Trou de serviette en forme de U
	Agujero en forma de U.
	U型ホール

Extremity Drape



- Ref no.:1.9.1.8
- Material:45-50gsm SMS, SMMS
- With circular hole retractable film around incise area
- Size : 80*100CM
- Packing: 1PC/POUCH , 100PCS/CTN

	Extremität Drapierung
	La serviette de trou
	Toalla de manos y pies
	両手足の穴

Laparotomy Drape




- Ref no.: 1.9.1.9
- Material:45gsm SMMS with fenestration incise film
- With adhesive tape around the fenestration
- Size : 80*100cm
- Packing: 1 set/pouch , 6 sets/ctn

	Laparotomie Drapierung
	Trou de serviettes abdominales
	Toalla abdominal
	腹部穴

C-section Drape



- Ref no.: 1.9.110
- Material: 45gsm SMMS, blue, with fluid collecting pouch, incise film
- With triangular/round fenestrated incise area and circular fluid collection pouch and filtration sponge
- Size: 260*300cm
- Packing: 1 pcs/pouch, 20 pouchs/ctn

	C-Abschnitt Drapierung
	La césarienne trou de serviette
	La cesárea toalla
	帝王切開洞

Angiography Drape






- Ref no.: 1.9.111
- Material: glued with surgery film, water-absorbent material or outside membrane, no woven fabric,
- With Two holes or three holes
- Size:customized
- Packing: 1pcs/pouch, 100pcs/ carton or customized

	Angiografie Drapierung
	Angiographie de trou de serviette
	Angiografía toalla
	血管造影

Cardiovascular Drape



- Ref no.: 1.9.112
- Material: SMS medical angiography drape
- With full chest and femoral incise area, adjustable tube holder and instrument pouches
- Size:customized
- Packing: 1pcs/pouch, 100pcs/ carton or customized

	Kardiovaskulär Drapierung
	Trou de tissu cardio-vasculaire
	Toalla cardiovascular
	心血管心臟心臟洞

Half/Large and Utility Sheets



- Ref no.: 1.9.113
- Dapes sheets are effective add-ons to other drapes to help create a sterile field and complete your procedural draping needs
- Size: customized
- Packing: 1pcs/pouch, 100pcs/ carton or customized

	Üblich Blätter
	Trou de serviette universel
	En general la toalla
	ユニバーサルホール

Disposable Mayo Stand Cover



- Ref no.: 1.9.114
- Material: PE film+PP nonwoven
- Size: 78cmx145cm, 80cmx140cm or Customized
- Packing: customized




	Einmal Mayo Stand Abdeckung
	Mayo unique
	Mayonesa simple
	梅奥羽書

Table Cover




- Ref no.: 1.9.115
- Material: PE or SMS with or without reinforcement
- High quality
- Size: 150*190CM or customized
- Packing: 30PCS/CTN
- Non sterile

	Wallpapers
	Fonds d'écran
	Fondos de pantalla
	壁紙


Disposable Surgical Urology Drape







- Ref no.: 1.9.116
- Material:45-55gsm SMS or SMMS
- High quality
- Size: 150*240CM or customized
- Packing: 30PCS/CTN

	Einmal chirurgisch Urologie Drapierung
	Serviette de trou d'urologie
	La Sección de la toalla de la cueva de urología
	泌尿科洞


Utility Drape



- Ref no.: 1.9.117
- Material: Impervious PE +absorbent tissue paper or Impervious PE +Non woven material
- EO/Radiation
- Machine folding with fan
- Color: blue or green
- Size: customized
- Packing: customized

	Üblich Drapierung.
	General de deux couches de trou de serviette
	El agujero de la capa de dos toallas
	G M二段穴


Fenestrated Drape



- Ref no.: 1.9.118
- Material: Impervious PE+absorbent tissue paper, with Hole: oval, roundness or rectangle
- Impervious PE+Non woven
- EO/Radiation
- Size: customized
- Packing: customized

	Fenestrated Drapierung
	Fenestrated Drapé
	Fenestrado cubrir
	開穴巾

Utility Drape with Adhesive and Hole



- Ref no.: 1.9.119
- Material: Impervious PE +Viscose, hole:8*8cm or Impervious PE+absorbent
- PP Diameter:10cm
- With Medical double-sided tape
- Sterile
- Size:50*50cm or customized
- Packing: customized

	Üblich Drapierung mit Klebstoff und fenester
	General avec de la colle de la serviette.
	Abrir agujeros con pegamento universal.
	汎用ゴム開穴巾



OFFICE:

**Jumeirah Lake Towers, Cluster G, JBC 1
Office 3801, Dubai, United Arab Emirates**

**info@gsticmedical.com
www.Qcaremd.com**