

2022




SURGICAL DRAPES AND PACKS

Surgical Basic Universal Drape Set



- Ref no.:1.9.2.1
1. Towel, 30x40cm, 4pcs
 2. Reinforced Gown, L, 2pcs
 3. Mayo Table Cover, 80x150cm, 1pc
 4. Adhesive Drape, 75x100cm, 2pcs
 5. Op-Tape, 10x50cm, 2pc

| | |
|---|--|
|  | Chirurgische allgemeine Drapierung Set |
|  | Ensemble de chirurgie générale |
|  | La cirugía general |
|  | 手術汎用セット |

Arthroscopy Pack Set



- Ref no.:1.9.2.2
1. BackTableCover, 150*200cm
 2. Mayo Stand Cover, 140*80cm
 3. Drape with tape, 200*15cm
 4. Plastic probe cover, 246*13cm
 5. Surgical stockinet, 22*115cm
 6. Plastic cover, 40*70cm
 7. Hand towels *2, 80*38cm
 8. OP tape*2, 50*10cm
 9. Kneearthroscopy drape, 330*200cm
 - 10.Surgical gown *3

| | |
|---|---------------------------------------|
|  | Arthroskopie Pack Set |
|  | Arthroscopie Pack jeu |
|  | Conjunto de paquete de la artroscopia |
|  | 関節鏡パックセット |

Angiography Drape Set



- Ref no.: 1.9.2.3
- PE film+PP nonwoven+airlaid paper
 - Two holes or three holes
 - Customized Accessories

| Name: | Size: | Quantity: |
|--------------------|-----------|-----------|
| Towel | 30*40cm | 4pcs |
| Reinforced Gown | L | 2pcs |
| Mayo Table Cover | 80*150cm | 1pc |
| Adhesive Drape | 75*100cm | 2pcs |
| OP-Tape | 10*50cm | 2pcs |
| Sterilization Wrap | 100*100cm | 1pc |

| | |
|---|------------------------------|
|  | Angiografie Drapierung Set |
|  | Angiographie drapé ensemble |
|  | Angiografía conjunto cortina |
|  | 血管造影用ドレープセット |

Laparotomy Pack Set



- Ref no.: 1.9.2.4
- 1PC BACK TABLE COVER 150*240CM
- 5PCS SURGICAL GOWN "L"
- 1 PC MAYO STAND COVER 60*140CM
- 4PCS UTILITY DRAPE WITH ADHESIVE EDGE 40*65CM
- 1PC LAPAROTOMY DRAPE 190*320CM
- 4PCS HAND TOWEL





| | |
|---|---------------------------------|
|  | Laparotomie Pack Set |
|  | Set pack laparotomie |
|  | Conjunto paquete de laparotomía |
|  | 開腹パックスセット |

C-section Drape Set



- Ref no.: 1.9.2.5

| Name: | Size: | Quantity: |
|--------------------|-----------|-----------|
| Towel | 30*40cm | 4pcs |
| Reinforced Gown | L | 2pcs |
| Mayo Table Cover | 80*150cm | 1pc |
| Adhesive Drape | 75*100cm | 2pcs |
| OP-Tape | 10*50cm | 2pcs |
| Sterilization Wrap | 100*100cm | 1pc |

| | |
|---|---------------------------------|
|  | C-Abschnitt Drapierung Set |
|  | C-section jeux de drapé |
|  | C-sección conjuntos de drapeado |
|  | 帝王切開セット |

Disposable TUR Surgical Set



- Refno.:1.9.2.6
- Specifications
 - Table cover 150x200cm Laminate
 - T.U.R. Drape 150/240x115cm SMS/Spunlace/Laminate
 - Legging 60x110cm SMS/Laminate
 - O.P. Tapes 10x50cm SMS/Spunlace
 - Hand towels 40x40cm paper

| | |
|---|------------------------------------|
|  | Einmal TUR chirurgisch Set |
|  | Set chirurgicale jetable TUR |
|  | Conjunto quirúrgica desechable TUR |
|  | 使い捨てTUR 手術用セット |

Disposable Cesarean Surgical set



- Ref no.:1.9.2.7

| Name: | Size: | Quantity: |
|--------------------|-----------|-----------|
| Towel | 30*40cm | 4pcs |
| Reinforced Gown | L | 2pcs |
| Mayo Table Cover | 80*150cm | 1pc |
| Adhesive Drape | 75*100cm | 2pcs |
| OP-Tape | 10*50cm | 2pcs |
| Sterilization Wrap | 100*100cm | 1pc |

| | |
|---|---|
|  | Einmal Cesarean chirurgisch Set |
|  | Jetable ensemble chirurgical Césarienne |
|  | Conjunto quirúrgico desechable cesárea |
|  | 使い捨て帝王切開手術用セット |

Surgery Urology Drape Pack



- Ref no.: 1.9.2.8

| Name: | Size: | Quantity: |
|--------------------|-----------|-----------|
| Towel | 30*40cm | 4pcs |
| Reinforced Gown | L | 2pcs |
| Mayo Table Cover | 80*150cm | 1pc |
| Adhesive Drape | 75*100cm | 2pcs |
| OP-Tape | 10*50cm | 2pcs |
| Sterilization Wrap | 100*100cm | 1pc |

| | |
|---|-------------------------------------|
|  | Chirurgie Urologie Drapierung Pack |
|  | Urologie pack drap |
|  | Cirugía paquete de cortina Urología |
|  | 手術泌尿器ドレーブパック |

Instrument Pack



- Ref no.:1.9.2.9

| Name: | Size: | Quantity: |
|-----------------------|-----------|-----------|
| PE instrument housing | 26x12.5cm | 2pcs |
| SMS sleeve | 210 x32cm | 1pc |
| Sterile bag | 35x30cm | 1pc |

| | |
|---|----------------------------|
|  | Instrumentenabdeckung |
|  | Couverture de l'instrument |
|  | Cubierta del instrumento |
|  | 計器カバー |

Universal Pack




- Ref no.:1.9.2.10

| Name: | Size: | Quantity: |
|------------------|-----------|-----------|
| Mayo stand cover | 80x145cm | 1pc |
| Side drape | 75x90 | 2pcs |
| Top drape | 150*240cm | 1pc |
| Bottom drape | 150*200cm | 1pc |
| OP-Tape | 10*50cm | 1pc |
| Hand towels | 30*40cm | 4pcs |
| Table cover | 160*190cm | 1pc |

| | |
|---|----------------------------|
|  | Allgemeine Chirurgie-Kits |
|  | Kits chirurgicaux généraux |
|  | Kits de cirugía general |
|  | 一般的な手術用キット |


Extremity Pack



- Ref no.:1.9.2.11
- Suture bag, 1pc
- Astound gown, impervious reinforced, X-large, 1pc
- Absorbent towel, 15 1/2 in.x22 in. 1pc
- Mayo stand cover, reinforced poly, 23 in x 55in 1pc
- Plastic sheet, fan-folded, 44 x 60 in. 1pc
- Drape towel, adhesive, nonabsorbent 3pcs
- Outer wrap 50 x 90 in.

| | |
|---|---------------------------|
|  | Allgemeine Chirurgie-Kits |
|  | Pack extrémité |
|  | Paquete de extremidad |
|  | 四肢パック |


Head and Neck Packs







- Ref no.:1.9.2.12
- 1pc head/bar drape, adhesive
- 1pc split sheet drape, reinforced, adhesive
- 1pc standard scrub nurse gown, large with towel in outer wrap
- 1pc standard gown, sterile back, large
- 2pcs Drape towel, adhesive, nonabsorbent
- 1pc Absorbent towel
- 1pc Mayo stand cover, reinforced poly, 23
- 1pc Suture bag
- 1pc Outer wrap, reinforced poly table cover, 44 x 90 in.

| | |
|---|--------------------------|
|  | Kopf-und Halspackungen |
|  | Packs de tête et du cou |
|  | Packs de cabeza y cuello |
|  | 頭頸部パック |

Cardio Vascular Set Sterile




- Ref no.:1.9.2.13
- 1pc Cardio-vascular drape 3-anti SMS 50gsm
- 1pc Table cover 150*200cm, 45gsm PP+ 60gsm PE
- 1pc Mayo cover, 80*140cm: PE, reinforcement: 140*60cm
- 1pc Drape sheet 150*175cm
- 2pc OP tape 10*50cm
- 4pc Towel 33*33cm

| | |
|---|----------------------------|
|  | Herz-Kreislauf -Set |
|  | Ensemble cardio-vasculaire |
|  | Conjunto cardiovasculares |
|  | 心血管セット |

Craniotomy Pack



- Ref no.:1.9.2.14
- Absorbent 35X22in
- Clear anesthesia panels
- Fenestration,oval with full incise 11x7in
- Fluid control drawstring pouch with screen and port
- Velcro cord-holding tabs

| | |
|---|----------------------------|
|  | Herz-Kreislauf -Set |
|  | Ensemble cardio-vasculaire |
|  | Conjunto cardiovasculares |
|  | 心血管セット |

Ophthalmic Pack-sterile



- Ref no.:1.9.2.15
- 1pc x hand towel, 30*40cm,
- 1pc x table cover,25gsm PE 150*200cm,reinforcement: 22gsm hydrophilic PP,80*150cm
- 1pc x Ophthalmic drape, hydrophilic PP+ PE100*120cm, aperture 8cm with liquid collection pouch 22*20CM & aluminum strip

| | |
|---|------------------|
|  | Augenpackung |
|  | Pack ophtalmique |
|  | Pack oftálmica |
|  | 眼科パック |



OFFICE:

**Jumeirah Lake Towers, Cluster G, JBC 1
Office 3801, Dubai, United Arab Emirates**

**info@gsticmedical.com
www.Qcaremd.com**