

2022



SURGICAL DRAPES AND PACKS

IV Start Kit



- Ref no.:1.9.3.1
- 1 x Alcohol pad
- 1 x Iodine pad
- 2x Nonwoven swabs.2"x2", 4-ply
- 1x Transparent Dressing
- 1x PE Tape roll 3/4"*18"
- 1x Tourniquet,Latex-free. 3/4"*18"
- 1 Pair Nitrile gloves
- 1x Patient ID Label

	Wundverband -Kits
	Kits de pansement
	Apósitos para heridas
	創傷被覆材パッケージ

Basic IV Start Kit



- Ref no.: 1.9.3.2
- 1x Alcohol pad
- 1x chlorhexidine pad
- 1x iodine pad
- 2x Nonwoven swabs.2"*2", 4-ply
- 1x Transparent Dressing
- 1x PE Tape roll 3/4"*18"
- 1x Tourniquet,Latex-free
- 1x Patient ID Label

	Allgemeine Wundverband -Kits
	Kits de pansement générales
	Apósitos para heridas generales
	ドレッシングパック

Basic Dressing Pack-Sterile



- Ref no.: 1.9.3.3
- 1x Premium Dressing Tray
- 6x Non Woven Swabs
- 1x Forceps Yellow
- 2x Forceps Green
- 1x Towel
- 1x PE wrapper
- 1x PE+paper bag
- Packing:1set/ PE+paper pouch, 100sets/ctn

	Allgemeine Dressing-Kits
	Kits générales de chambre
	Kits de vestidores generales
	一般的なドレッシングキット

Standrad Dressing Pack-Sterile



- Ref no.: 1.9.3.4
- 1x Premium Dressing Tray
- 1x Forceps Yellow
- 6x Non Woven Swabs
- 1x Towel Blue Polylined
- 2x Forceps Green
- 1x White Plastic Field

- Packing:1set/ PE+paper pouch,
60sets/ctn

	Standard-Dressing-Kits
	Assortiments de pansements standard
	Apósitos estándar
	標準ドレッシングキット

Dressing Pack-Sterile



- Ref no.: 1.9.3.5
- 1x Blue wrapper 49.7*78cm 20gsm film
- 1x Drape 49.3 x 49cm, 21gsm paper+14gsm PE
- 1x Compartmented tray 14x12x2cm
- 5x non-woven swab 7.5x7.5cm x 4ply 30gsm spunlace
- 8x cotton wool ball 0.5gsm
- 2x green forceps
- 1x yellow forcep
- Packing: 1set/ PE+paper pouch,
100sets/ctn

	Dressing Pack
	Kits de pansements
	Apósitos paquete
	ドレッシングキット

Dressing Pack-Sterile



- Ref no.: 1.9.3.6
- 1x blue wrapper 45*45cm
- 1x Tray 7.2*10.2*3.2cm
- 4x Forceps Blue
- 5x non-woven roll
- 1x yellow disposable bag
- 1x Towel,55*55cm
- 5x nonwoven swabs, 4-ply,7.5*7.5cm
- Packing:1set/ PE+paper pouch,
100sets/ctn

	Dressing Pack
	Kits de pansements
	Apósitos paquete
	ドレッシングキット

Skin Excision Kit



- Ref no.: 1.9.3.7
1. Paper Drape 80*50cm, 1 pc
 2. Scalpel with Blade #11, 1pc
 3. Non-Woven Swabs 4-ply 7.5*7.5cm, Spunlace 30gsm, 2pcs
 4. Non-Woven Balls 0.5gsm, 5pcs
 5. Yellow Forceps 10.8cm, pc
 6. Blue Forceps 10.8cm, pc
 7. Field PE Film 50*50cm, 1pc
 8. Packing: EO sterilization, 1kit/bag, 100bag/ctn, 66*30*33cm

	Hautexzision KIT
	EXCISION DE LA PEAU KIT
	KIT escisión de la piel
	皮膚切除KIT

Staple Removal Kit



- Ref no.: 1.9.3.8
1. Staple removal 10cm, stainless, 1 pc
 2. Gauze Sponge 40's 7.5*7.5cm-8ply, 1pc
 3. Packing: 1kit/blister card, 100kits/ctn 53*34*27cm

	HEFT Entfernungskit
	STAPLE RETRAITS KIT
	Kit de eliminación de ENGRAPE
	STAPLE 除去キット

Suture Removal Kit



- Ref no.: 1.9.3.9
1. Littauer Scissors 4.5", 1 pc
 2. Metal Forceps 4", 1pc
 3. Gauze Sponge, 1pc
 4. Packing: 1kit/blister card, 50kits/ctn 28*28*17cm

	SUTURE Entfernungskit
	SUTURE RETRAITS KIT
	Estuche de sutura EXTRACCIÓN
	縫合糸除去キット

Male Circumcision Kit



- Ref no.: 1.9.3.10

1. Multipurpose container tray- 26cmx18cmx3cm, 1 pc
2. O-Drape, hole diameter 5cm, 100cmX75cm, 1pc
3. Gauze swab, 10cmx10cm, 12ply, 20pcs
4. Paraffin gauze 10cmx10cm, 1pc
5. Syringes, 10mls, 1pc
6. Injection needle 21G and 23%, 2pcs
7. Sutures, braded/absorbable w/needle, 3/0, 75cm---2pcs
8. Surgical gloves, size size 7.5 and 8, 2 pairs
9. Isopropylalcohol swab 3cmx6cm, 2 pcs
10. Plastic Apron, 80cmx130cm, 2pcs
11. Surgical paper tape 1.2cmx1cm, 1pc
12. Examination gloves, Large, 1 pairs
13. Needle holder 14cm, 1pc
14. Suture scissor 14cm, 1pc
15. Non toothed forceps, Steel, 13cm, 1pc
16. Mosquito clamp, curved,14cm, 1pc
17. Mosquito clamp, straight, 14cm, 1pc
18. Circumcision forceps, haemostatic, cross-clamp
19. Scalpel handle, retractable, no 23m 1pc
20. Surgical crepe paper 60cx60cm, 1pc
21. EO indicators, 1pc
22. Packing: 1kit/box, 20boxes/ctn, 45*30*40cm.

	MALE BESCHNEIDUNG KIT
	CIRCONCISION MASCULINE KIT
	KIT DE LA CIRCUNCISIÓN MASCULINA
	男性の割礼キット

Male Circumcision Kit



- Ref no.:1.9.3.11

1. Multipurpose container tray- 26cmx18cmx3cm, 1 pc
2. O-Drape, hole diameter 5cm, 100cmX75cm, 1pc
3. Gauze swab, 10cmx10cm, 12ply, 20pcs
4. Paraffin gauze 10cmx10cm, 1pc
5. Syringes, 10mls, 1pc
6. Injection needle 21G and 23%, 2pcs
7. Sutures, braded/absorbable w/needle, 3/0, 75cm---2pcs
8. Surgical gloves, size size 7.5 and 8, 2 pairs
9. Isopropylalcohol swab 3cmx6cm, 2 pcs
10. Plastic Apron, 80cmx130cm, 2pcs
11. Surgical paper tape 1.2cmx1cm, 1pc
12. Examination gloves, Large, 2 pairs
13. Needle holder 14cm, 1pc
14. Suture scissor 14cm, 1pc
15. Non toothed forceps, Steel, 13cm, 1pc
16. Metzenbum Scissors,curved 15cm, 1pc
17. Metzenbum Scissors,straight, 15cm, 1pc
18. Circumcision forceps, haemostatic, cross-clamp 20cm, 1pc
19. Scalpel handle, retractable, no 23m 1pc
20. Surgical crepe paper 60cx60cm, 1pc
21. EO indicators, 1pc
22. Packing: 1kit/box, 20boxes/ctn, 45*30*40cm.

	MALE BESCHNEIDUNG KIT
	CIRCUNCISION MASCULINE KIT
	KIT DE LA CIRCUNCISIÓN MASCULINA
	男性の割礼キット



OFFICE:

**Jumeirah Lake Towers, Cluster G, JBC 1
Office 3801, Dubai, United Arab Emirates**

**info@gsticmedical.com
www.Qcaremd.com**