

2022



**PROTECTIVE
DISPOSABLES**

Isolation Gown with Elastic Cuff



- Ref no.: 1.1.3.4.1
- Material: PP Nonwoven or SMS and so on
- With elastic cuff neck and 2 waist tape
- Size: 110*120cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	Isolation Kleid
	Robe d'isolement
	Bata de aislamiento
	アイソレーションガウン

Isolation Gown with Knitted cuff



- Ref no.: 1.1.3.4.2
- Material: PP Nonwoven or SMS and so on
- With knitted cuff neck and 2 waist tape
- Size: 110*120cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	Isolation Kleid
	Robe d'isolement
	Bata de aislamiento
	アイソレーションガウン

Surgical Gown



- Ref no.: 1.1.3.7.1
- Material: SMS
- With knitted cuff Velcro neck and 2 waist tape
- Size: 110*120cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	OP-Kittel
	Blouse chirurgicale
	Bata quirúrgica
	手術衣

Reinforced Isolation Gown



- Ref no.: 1.1.3.4.3
- Material: PP+PE
- With knitted cuff, Velcro neck and waist tape.
- Reinforced with PP+PE for the whole gown
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	OP-Kittel verstärkt
	Blouse chirurgicale renforcée
	Bata quirúrgica reforzada
	手術用ガウンを強化します

Reinforced Isolation Gown



- Ref no.: 1.1.3.4.4
- Material: PP+PE
- With elastic cuff, Velcro neck and waist tape.
- Reinforced with PP+PE for sleeve and front area.
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	OP-Kittel verstärkt
	Blouse chirurgicale renforcée
	Bata quirúrgica reforzada
	手術用ガウンを強化します

Reinforced Surgical Gown



- Ref no.: 1.1.3.7.3
- Material: SMS
- With knitted cuff, Velcro neck and 4 waist tape.
- Reinforced with PP+PE in front of gown
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	OP-Kittel verstärkt
	Blouse chirurgicale renforcée
	Bata quirúrgica reforzada
	手術用ガウンを強化します

Reinforced Isolation Gown by Ultrasonic



- Ref no.:1.1.3.4.3.1
- Material: PP+PE
- With elastic cuff, neck and waist tape. by ultrasonic.
- Reinforced with PP+PE for the whole gown
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	OP-Kittel verstärkt
	Blouse chirurgicale renforcée
	Bata quirúrgica reforzada
	手術用ガウンを強化します

Surgical Gown by Ultrasonic



- Ref no.:1.1.3.7.1.1
- Material: SMS
- With knitted cuff Velcro neck and 2 waist tape, by ultrasonic.
- Size: 110*120cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	OP-Kittel verstärkt
	Blouse chirurgicale renforcée
	Bata quirúrgica reforzada
	手術用ガウンを強化します

Reinforced Surgical Gown by Ultrasonic



- Ref no.:1.1.3.7.3.1
- Material: SMS
- With knitted cuff, Velcro neck and 4 waist tape, by ultrasonic.
- Reinforced with PP+PE in front of gown
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	OP-Kittel verstärkt
	Blouse chirurgicale renforcée
	Bata quirúrgica reforzada
	手術用ガウンを強化します

Patient Gown



- Ref no.: 1.1.3.4.1.2
- Material: PP Nonwoven or SMS and so on
- without sleeve neck and 2 waist tape
- Size: 120*140cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	Patientenhemd
	Robe du patient
	Bata de paciente
	患者ガウン

Childbirth Shirt



- Ref no.: 1.1.3.10
- Material: PP Nonwoven
- Size: 180*110cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	Child Schönheit rth auf rt
	La beauté de rth de Chi à la température ambiante
	Child rth belleza a rt
	出産のシャツ

Patient Exam Capes



- Ref no.: 1.1.3.10.1
- Material: PP Nonwoven
- No sleeve, neck and waist with tie
- Size: 75*52cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

	Prüfung Umhänge
	Capes d'examen
	Cabos del examen
	試験の岬

Spunlace Surgical Gown



- Ref no.:1.1.3.7.4
- Material: Spun lace wood pulp
- With knitted cuff, Velcro neck and 4 waist tape
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: blue or all color can be customized

	Spunlaced-OP-Kittel
	Dentelle filée blouse chirurgicale
	Bata quirúrgica de encaje hilado
	スパンレース手術用ガウン

Scrub Suits



- Ref no.:1.1.3.8
- Material: PP Nonwoven or SMS
- Short sleeve With V neck or round neck with or without pocket, Elastic or tie waist
- Top Size: 70*110cm pant size: 105*100cm, size can be customized
- Packing: 1pcs/PE bag, 50bags/ctn, can be customized
- Color option: blue or all color can be customized

	Bereichskleidung
	Costumes gommage
	Trajes de matorral
	スクラブスーツ

Exam Gown



- Ref no.:1.1.3.9
- Material: PP Nonwoven or SMS
- Short sleeve or no sleeve with or without pocket,
- Top Size: 110*120cm, size can be customized
- Packing: 1pcs/PE bag, 50bags/ctn, can be customized
- Color option: blue or all color can be customized

	Prüfung Kleid
	Robe examen
	Bata de examen
	試験のガウン

Coverall



- Ref no.:1.1.3.3
- Material: PP Non woven or SMS
- With hood or without hood, zipper, elastic cuff and bottom, waist
- Size: 120*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or black color can be customized

	Overall
	Coverall
	Sobretudo
	カバーオール

Coverall Type 5/6



- Ref no.:1.1.3.3.3
- Material: PP Non woven + Breathable film
- With hood, zipper, elastic cuff and bottom, waist
- Size: 120*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or blue, color can be customized

	Overall
	Coverall
	Sobretudo
	カバーオール

Coverall Type 4/5/6



- Ref no.:1.1.3.3.4
- Material: PP Non woven + Breathable film
- Elastic in waist area providing body-fitting of a coverall, elastic in cuffs, legs and hood, loose shape with gore in groin area, three-part hood, self-adhesive cover of a zipper, adhesive chin cover, with elastic thumb loops
- Size: 120*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or blue, color can be customize

	Overall
	Coverall
	Sobretudo
	カバーオール

Lab Coat



- Ref no.:1.1.3.5
- Material: PP Non woven or SMS
- With collar, Elastic or knitted cuff Velcro or snap closure
- Size: 120*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or black color can be customized

	Laborkittel
	Blouse de laboratoire
	Bata de laboratorio
	白衣

Lab Coat



- Ref no.:1.1.3.5.1
- Material: PP Non woven or SMS
- With Knitted collar, Elastic or knitted cuff Velcro or snap closure
- Size: 120*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or black color can be customized

	Laborkittel
	Blouse de laboratoire
	Bata de laboratorio
	白衣

Non Woven Kimono



- Ref no.:1.1.3.14
- Material: PP Non woven
- With one tie With one pocket in short sleeve
- Size: 100*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or black color can be customized

	Non-woven Kimono
	Kimono non tissé
	Kimono no tejidas
	不織布着物

Non Woven Shirt



- Ref no.:1.1.3.11
- Material: PP Non woven/ PP+PE
- Short sleeve and v neck
- Size: 50*120cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: White, color can be customized

	Die do Oh p Theorie Ven
	Ven le do de la théorie de p Oh
	Ven el do de Teoría p Oh
	不織布エプロン

Non Woven Pant



- Ref no.:1.1.3.12
- Material: PP Non woven
- With elastic or tie waist Elastic bottom
- Size: 105*110cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or black color can be customized

	Non-woven Verräter
	Non tissé traître
	No tejida traidor
	不織布裏切り者

Non Woven Apron



- Ref no.:1.1.3.16
- Material: PP Non woven/ PP+PE
- With tie in waist
- Size: 75*100cm size can be customized
- Packing: 50pcs/box, 8boxes/ctn, can be customized
- Color option: White, color can be customized

	Die do Oh p Theorie Ven
	Ven le do de la théorie de p Oh
	Ven el do de Teoría p Oh
	不織布エプロン

Non Woven Underwear



- Ref no.:1.1.3.1.1
- Material: PP Non woven
- With elastic
- Size: 52*25cm size can be customized
- Packing: 1pc/PE bag, 1000bags/ctn, can be customized
- Color option: white or black color can be customized

	r r Hardware in Ja Ven die Do von
	r r Hardware dans Ouais Ven le do de
	La ropa interior no tejida
	不織布下着

Non Woven Underwear



- Ref no.:1.1.3.1.2
- Material: PP Non woven
- With elastic
- Size: 44*16cm size can be customized
- Packing: 1pc/PE bag, 1000bags/ctn, can be customized
- Color option: white or black color can be customized

	r r Hardware in Ja Ven die Do von
	r r Hardware dans Ouais Ven le do de
	La ropa interior no tejida
	不織布下着

Non Woven Bra



- Ref no.:1.1.3.15
- Material: PP Non woven
- With elastic
- Size: 135*20cm, size can be customized
- Packing: 50pc/PE bag, 500bags/ctn, can be customized
- Color option: white or black color can be customized

	Non-Woven-BH
	Non tissé soutien-gorge
	Sujetador no tejida
	不織布ブラジャー

PE Sleeves



- Ref no: 1.1.3.17
- Made of PE. With comfortable elastic opening
- Waterproof and dustproof
- Light weight, convenience, flexible, economical, soft
- Available in different thickness and color
- Size and packing can be customized

	PE-Mantel
	PE manchon
	Manga de PE
	PE スリーブ

PE Apron



- Ref no: 1.1.3.18
- Made of LDPE, HDPE
- Size: 80*120cm, 80*125cm
- Color: white
- Packing: 1000pcs/ctn

	PE Schürze
	PE Tablier
	PE delantal
	PE エプロン

PE Visit Gown



- Ref no: 1.1.3.19
- Made of PE.
- Size: 120*145cm
- Color: white
- Packing: 100pcs/ctn

	PE Schürze
	PE Tablier
	PE delantal
	PE エプロン

CPE Gown



- Ref no: 1.1.3.20
- Made of CPE.
- Size:120*198cm or other size
- Color: blue
- Packing:100pcs/ctn

	CPE Schürze
	CPE Tablier
	CPE delantal
	CPE エプロン

Dental Bib



- Ref no: 1.1.3.21
- Size: 33*45cm
- Wood pulp paper + water-proof PE backing
- 2 plies or 3plies available
- Packing: 125pcs/ bag 4bags/ctn
- Color option: different color

	Dental bib
	Babero dental
	Den era l belleza b
	歯科よだれかけ

Dental Apron Finland Type



- Ref no: 1.1.3.22.1
- Made of wet strong paper and water proof PE backing
- With a smart strip to guarantee comfort wear
- with pocket in front
- Packing: 300pcs/ctn
- Color option: blue white

	Dental Schürze Finnland Art
	Type Tablier dentaire Finlande
	Tipo Delantal Dental Finlandia
	歯科エプロンフィンランド式

Dental Apron Finland Type



- Ref no: 1.1.3.22.2
- Made of wet strong paper and water proof PE backing
- With adhesive sticker to guarantee comfort wear
- with pocket in front
- Packing: 300pcs/ctn
- Color option: blue white

	Dental Schürze Finnland Art
	Type Tablier dentaire Finlande
	Tipo Delantal Dental Finlandia
	歯科エプロンフィンランド式

Dental Apron Italy Type



- Ref no: 1.1.3.23.1
- Made of wet strong paper and water proof PE backing - With adhesive sticker to guarantee comfort wear
- With overhead
- Packing: 300pcs/ctn
- Color option: blue white pink

	Dental Schürze Italien Art
	Type Tablier dentaire Italie
	Tipo Delantal Dental Italia
	歯科エプロンイタリアタイプ

Dental Apron Italy Type



- Ref no: 1.1.3.23.2
- Made of wet strong paper and water proof PE backing
- With adhesive sticker to guarantee comfort wear
- With pocket backside
- Packing: 300pcs/ctn
- Color option: blue white pink

	Dental Schürze Italien Art
	Type Tablier dentaire Italie
	Tipo Delantal Dental Italia
	歯科エプロンイタリアタイプ

Dental Apron German Type



- Ref no: 1.1.3.24.1
- Made of wet strong paper and water proof PE backing
- With a smart strip to guarantee comfort wear
- With pocket backside
- Packing: 300pcs/ctn
- Color option: blue white

	Dental Schürze deutschen Typ
	Tablier dentaire de type allemand
	Delantal Dental tipo alemán
	歯科エプロンドイツタイプ

Dental Apron German Type



- Ref no: 1.1.3.24.2
- Made of wet strong paper and water proof PE backing
- With a smart strip to guarantee comfort wear
- With pocket backside
- Packing: 300pcs/ctn
- Color option: blue white

	Dental Schürze deutschen Typ
	Tablier dentaire de type allemand
	Delantal Dental tipo alemán
	歯科エプロンドイツタイプ

Dental Apron American Type



- Ref no: 1.1.3.25.1
- Made of wet strong paper and water proof PE backing
- With a smart strip to guarantee comfort wear
- With tie, no pocket
- Packing: 300pcs/ctn
- Color option: blue white

	Dental chürze amerikanische Art
	Type dentaire Tablier américain
	Tipo Delantal Dental Americana
	歯科エプロンアメリカンタイプ

Dental Apron American Type



- Ref no: 1.1.3.25.2
- Made of wet strong paper and water proof PE backing
- Overhead, no pocket
- Packing: 300pcs/ctn
- Color option: blue white

	Dental chürze amerikanische Art
	Type dentaire Tablier américain
	Tipo Delantal Dental Americana
	歯科エプロンアメリカンタイプ

Dental Apron Type



- Ref no: 1.1.3.26.1
- Made of wet strong paper and water proof PE backing
- Overhead, no pocket
- Packing: 300pcs/ctn
- Color option: blue white

	Dental Schürze Art
	Type Tablier dentaire
	Tipo Delantal Dental
	歯科エプロンイスラエルタイプ

Dental Apron Type



- Ref no: 1.1.3.26.2
- Made of wet strong paper and water proof PE backing
- overhead, no pocket
- Packing: 300pcs/ctn
- Color option: blue white

	Dental Schürze Art
	Type Tablier dentaire
	Tipo Delantal Dental
	歯科エプロンイスラエルタイプ

Dental Apron Spain Type



- Ref no: 1.1.3.27.1
- Made of wet strong paper and water proof PE backing
- There are two tapes on the shoulder. It's convenient.
- Packing: 300pcs/ctn

	Spanien Lätzchen
	Espagne Bibs
	España delantal dental
	歯科スペインビブ

Special Dental Apron



- Ref no: 1.1.3.28.1
- Made of wet strong paper and water proof PE backing
- Size:30*40cm overhead, no pocket
- Packing: 300pcs/ctn
- Color option: white

	Dental Schürze
	Tablier dentaire
	Delantal dental
	歯科エプロン

Special Dental Apron



- Ref no: 1.1.3.28.2
- Made of wet strong paper and water proof PE backing
- Size:40*76cm,
- Packing: 300pcs/ctn
- Color option: white, blue

	Dental Schürze
	Tablier dentaire
	Delantal dental
	歯科エプロン

Special Dental Apron



- Ref no: 1.1.3.28.3
- Made of wet strong paper and water proof PE backing
- Size: 40*124cm,
- Packing: 300pcs/ctn
- Color option: Blue,white

	Dental Schürze
	Tablier dentaire
	Delantal dental
	歯科エプロン

Special Dental Apron



- Ref no: 1.1.3.28.4
- Made of wet strong paper and water proof PE backing
- Size: 37*63cm.pocket with paper outside
- Packing: 300pcs/ctn
- Color option: Blue,white

	Dental Schürze
	Tablier dentaire
	Delantal dental
	歯科エプロン

Special Dental Apron



- Ref no: 1.1.3.28.5
- Made of wet strong paper and water proof PE backing
- Size: 37*46,37*62cm.pocket with paper outside
- Packing: 300pcs/ctn
- Color option: Blue,white

	Dental Schürze
	Tablier dentaire
	Delantal dental
	歯科エプロン

VISITOR KIT



- Ref no.:1.1.3.29.1
- 1. Clip Cap: PP10gsm, Size 21" --1 pc
- 2. Shoe Cover: PP25gsm, 17*38cm, 2pcs
- 3. Face Mask: PP18/20/25gsm, 17.5*9.5cm --1pc
- 4. Isolation Gown: PP25gsm, 130*140cm --1pc
- Packing: 1kit/bag, 20bags/ctn, 29*25*30cm
- The components can be customized

	Besucher-Kit
	Kit de visiteur
	Kit de visitante
	訪問者キット



VISITOR KIT



- Ref no.:1.1.3.29.2
- 1. Clip Cap: PP10gsm, Size 21" --1 pc
- 2. Shoe Cover: PP25gsm, 17*38cm, 2pcs
- 3. Face Mask: PP18/20/25gsm, 17.5*9.5cm --1pc
- 4. Visitor gown: PP25gsm,130*140cm --1pc
- Packing: 1kit/bag, 20bags/ctn, 29*25*30cm
- The components can be customized

	Besucher-Kit
	Kit de visiteur
	Kit de visitante
	訪問者キット



VISITOR KIT



- Ref no.:1.1.3.29.3
- 1. Clip Cap: PP10gsm, Size 21" --1 pc
- 2. Shoe Cover: PP25gsm, 17*38cm, 2pcs
- 3. Face Mask: PP18/20/25gsm, 17.5*9.5cm --1pc
- 4. Visitor gown: PE30g/pc, 130*140cm --1pc
- Packing: 1kit/bag, 20bags/ctn, 29*25*30cm
- The components can be customized



	Besucher-Kit
	Kit de visiteur
	Kit de visitante
	訪問者キット

Mother KIT



- Ref no.:1.1.3.29.4
- 1. Clip Cap: PP10gsm, Size 21" --1 pc
- 2. Shoe Cover: PP25gsm, 17*38cm, 2pcs
- 3. Face Mask: PP18/20/25gsm, 17.5*9.5cm --1pc
- 4. Mother Shirt: Spunlace 50gsm+SMS 35gsm,130*140cm --1pc
- Packing: 1kit/bag, 20bags/ctn, 29*25*30cm
- The components can be customized



	Mutter-Kit
	Kit mère
	Kit de la madre
	母キット



OFFICE:

Jumeirah Lake Towers, Cluster G, JBC 1
Office 3801, Dubai, United Arab Emirates

info@gsticmedical.com
www.Qcaremd.com