

2022



# TAPES AND DRESSINGS

## Nonwoven Tape



- Ref no:1.2.1.1
- Material: microporous non-woven fabric and medical
- Good Permeable
- Low Sensitive
- Good adhesiveness to skin
- Size:1.25cm\*9.1m,2.5cm\*9.1m,5cm\*9.1m
- Packing can be customized

	Gewebtes Band
	Ruban tissé
	Cinta tejida
	織テープ

## PE Tape



- Ref no:1.2.1.2
- Made of PE material
- Porous, and breathable, transparent
- Moderate viscosity, low allergy
- Various sizes for choice
- Size: 1.25cm\*9.1m,2.5cm\*9.1m,5cm\*9.1m
- Packing can be customized

	PE-Band
	Bande PE
	Cinta de PE
	PE テープ

## Silk Tape



- Ref no:1.2.1.3
- Material:silk fabric
- Proper Viscidity,
- Good Permeable,
- Comfortable feel and no harm to skin
- Easily torn by hands
- Size: 1.25cm\*9.1m,2.5cm\*9.1m,5cm\*9.1m
- Packing can be customized

	Seidenband
	Ruban de soie
	Cinta de seda
	シルクテープ

## Heavy Elastic Sports Tape



- Ref no:1.2.1.4
- Material: cotton fabric coated with medical hot melt adhesive
- Sports support and protection and wound dressing
- Conforms easily to body contours and joints
- High tensile strength
- Size: 5cm\*4.5m,7.5cm\*4.5m,10cm\*4.5m
- Packing can be customized

	Sportband
	Bande de sport
	Cinta de los deportes
	スポーツテープ

## Double-sided Tape



- Ref no:1.2.1.5
- Material: Polyester
- Excellent temperature, water resistance
- Size: 1.25cm\*9.1m,2.5cm\*9.1m,5cm\*9.1m
- Packing can be customized

	Papier Klebeband
	Bande de papier
	Cinta de papel
	紙テープ

## Zinc Oxide Adhesive Plaster



- Ref no:1.2.1.6
- Material:cotton
- Natural rubber and zinc oxide adhesive
- Rasy to take
- Size: 5cm\*4.5m,7.5cm\*4.5m,10cm\*4.5m
- Packing can be customized

	Zinkoxid Band
	Bande d'oxyde de zinc
	Cinta de óxido de zinc
	酸化亜鉛テープ



## **OFFICE:**

**Jumeirah Lake Towers, Cluster G, JBC 1  
Office 3801, Dubai, United Arab Emirates**

**[info@gsticmedical.com](mailto:info@gsticmedical.com)  
[www.Qcaremd.com](http://www.Qcaremd.com)**